

Amvisc[®]

1.2% Sodium Hyaluronate

Vi
Viscoelastic

Amvisc[®] a pure sodium hyaluronate cohesive viscoelastic for all cataract procedures

Efficient Removal

Lasting Chamber Retention



Cohesive-Dispersive

Rheological properties

Composition	1.2% Sodium Hyaluronate
Buffered pH	ranging between 6.8 and 7.6
Molecular Weight	1-2.9 million Daltons
Viscosity	40,000 cP (at 1s ⁻¹ , 25°C)
Cohesion	high
Osmolality	320 mOsm

BAUSCH + LOMB

The benefits

- Ideal for the surgeon who performs cataract surgery using the planned extracaps technique.
- Its high viscosity provides optimal chamber maintenance during lens insertion and removal.
- Its high cohesive characteristic enables complete and efficient removal at the end of the procedure.

- Amvisc is available in 0.8ml single use pre-filled syringes
- Amvisc is a class IIb medical device (CE1275)
- Amvisc should be stored between 2°C to 8°C (36°F to 46°F) but allowed to reach room temperature (approx. 20 to 45 min) prior to use
- Protect from freezing.



Amvisc	0.8ml
Sterilisation	Sterile filtered, aseptically filled, secondary package sterile by EtO (Ethylene Oxide)
Cannula	27 Gauge
Order code	59081L

A complete family of Viscoelastic Solutions

Amvisc® Plus	The viscoelastic that adds the 'plus' of cohesive versatility for every stage of the cataract procedure. Amvisc Plus is molecularly engineered with a versatile range of cohesion that provides lasting chamber retention plus efficient removal at the end of the case. Cohesive versatility allows you to do what you want to do throughout the procedure, without the need for a second viscoelastic.	1.6% Sodium Hyaluronate	0.8 ml	60081L
OcuCoat	Is ideal for high volume anterior segment surgery. Formulated from low molecular weight 2% hydroxypropylmethylcellulose; allows rapid removal from the anterior chamber; has excellent coating properties making it ideal to use directly on the cornea to help protect delicate cells and maintain clarity throughout the procedure.	2% HPMC	1ml 1mlx6 2ml	CC050S CC065S CC100SL

Belgium/
Luxembourg Tel: +32 3 280 82 40
Fax: +32 3 280 82 59

Emerging Markets Tel: + 33 4 67 12 30 30
Fax: + 33 4 67 12 30 32

France Tel: +33 4 67 12 30 30
Fax: +33 4 67 12 30 31

Germany/
Switzerland/
Austria Tel: +49 30 33093 5431
Fax: +49 30 33093 5470

Italy Tel: +39 029 148 3851

Netherlands Tel: +31 20 6554555
Fax: +31 20 655 4640

Nordic Countries Tel: +46 8 616 95 00
Fax: +46 8 669 86 23

Portugal Tel: +351 214241425
Fax: +351 214241518

Spain Tel: +34 902 381 010
Fax: +34 902 250 310

South Africa Tel: +27 11372 5600

United Kingdom Tel: +44 20 8781 0000
Fax: +44 20 8781 0001

B+L



PREMED PHARMA KFT.
CÍM 2040 Budaörs, Gyár u. 2.
TELEFON 06 23 889 700
FAX 06 23 889 710
E-MAIL info@premedpharma.hu
WEB www.premedpharma.hu