

**VINNOA6**  
Simplicity Outside Intelligence Inside



**VINNO**

**VINNO Technology (Suzhou) Co., Ltd.**

5F, A Building, No.27 Xinfra Rd, Suzhou Industrial Park, 215123, China

Tel: +86 512 62873806

Fax: +86 512 62873801

Email: [vinno@vinno.com](mailto:vinno@vinno.com)

Website: [www.vinno.com](http://www.vinno.com)

VINNO reserves the rights to revise the technical specifications if needed.



**PREMEDI PHARMA KFT.**

**CÍM** 2040 Budaörs, Gyár u. 2.

**TELEFON** 06 23 889 700

**FAX** 06 23 889 710

**E-MAIL** [info@premedpharma.hu](mailto:info@premedpharma.hu)

**WEB** [www.premedpharma.hu](http://www.premedpharma.hu)

**VINNO**



**VINNO A6 is a brand-new Color Doppler Ultra-sound System. Adhering to VINNO's simplified design style and extremely smooth work flow, it will bring you a distinct diagnostic experience.**






# Simplified Innovation

## 1. Simplified Platform

- » Minimize a powerful platform into a compact size to meet your daily work.
- » Innovative cooling solution delivers comfortable and silent working environments.




Intergrate



Compact


## 2. Simplified UI



**Fast Setting Area**-Complete common setting under scanning status.  
**Outdoor and Indoor Solution**-Adjust monitor brightness according to the working environment.




**Quick Navigation Bar**-Browse long page fast by moving the trackball after click the interested area.  
**Back ICON**-Return to the top of the long page with only one click.




**Simple UI**-Select preset directly without reiterative page turning.  
**Magic Drag and Drop**-Adjust the sequence of measurement item easily by Magic Drag and Drop operation.

# Intelligent Workflow

## 1. Auto Image Optimization



off



on

## 2. Auto Calculation


### Auto OB

Just one click, the system will recognize the target automatically and the gestational age will be displayed there.



### Auto IMT

Just one click, the system will display max, average values of intima thickness and valid points at the same time.








## Excellent Image Quality

**VFusion:** Compounding technology together with RF-platform, VFusion helps to improve the image information, enhances contrast resolution and reduces the noise by scanning from different angles and combining those images.

**VSpeckle:** Reduces speckles in homogeneous tissue and enhance edge sharper than conventional image processing.

**VTissue:** Matches different acoustic velocity with different tissue automatically. It ensures the acquired image with more uniformity and accuracy.



# Beyond Your Expectation

- » 15.6 inch LED display
- » Monitor can be tilted up to 40 degrees
- » Detachable probe holder
- » Support up to 3 probe ports
- » battery lasting for about 1.5 hours
- » Support CWD, ECG, 4D function
- » Support WiFi/Bluetooth connectivity
- » Support transducer: A2-5C, A4-12L, A1-4P, A4-9E, A3-6D







## Enhance Your Skill

"Education is the passport to the future". Log on VINNO's Flyinsono APP and VCloud platform to study continuously and share your experiences.

